



ecoin.us

1-833-ECOINUS

1-833-326-4687

The eCoin Peripheral Neurostimulator (“eCoin”) is part of the eCoin Peripheral Neurostimulator System (“eCoin System”) including device accessories. Indications: The eCoin is intended to be used to treat urgency urinary incontinence in patients intolerant to or having an inadequate response to other more conservative treatments or who have undergone a successful trial of percutaneous tibial nerve stimulation. **Contraindications:** The eCoin is contraindicated for the following patients: • Poor Surgical Candidates: The eCoin should not be implanted in patients who are poor surgical candidates. Poor surgical candidates include those who have open wounds or sores on the lower leg or foot; had prior surgery in the implant area; had previous, unhealed trauma in the implant area; pitting edema ($\geq 2+$) in the lower leg; venous disease/insufficiency in the lower leg; arterial disease/insufficiency in the lower leg; vasculitis or dermatologic conditions in the lower leg, infections near the implantation site in the lower leg. • Patient cannot properly operate the Patient Controller Magnets and paper tape for use in the event of unintended or unwanted stimulation. **Warnings:** Full-body security scanners used by the TSA are considered safe in patients that have a stimulator. The eCoin is MRI Conditional. It is not safe to have a patient’s lower leg placed in an MRI machine. For a complete list of warnings, precautions, and potential adverse events, refer to the eCoin Physician or Patient Manual or visit www.ecoin.us/safety. Rx only.

eCoin® and Valencia Technologies are trademarks of Valencia Technologies Corporation, registered in U.S. and other countries. © 2022 Valencia Technologies. All rights reserved. 2221279834 Rev B